

Subject ID: _____

Date of Test: _____

Name/Initials of Tester: _____

Item	Item Name	(fill the section below that applies)							
		Right		Both	Left				
1	REACH FORWARD (reps)	<input type="checkbox"/>	Cannot do at all	Partial = hand at shoulder level, extended at least as far as knee; SEE MANUAL.	<input type="checkbox"/>	Cannot do at all			
		<input type="checkbox"/>	Partial completion of one		<input type="checkbox"/>	Partial completion of one			
			# reps/30 secs			# reps/30 secs			
2	REACH UP (reps)	<input type="checkbox"/>	Cannot do at all	*arm must be straight to count full repetition; Partial = hand above top of head; SEE MANUAL.	<input type="checkbox"/>	Cannot do at all			
		<input type="checkbox"/>	Partial completion of one		<input type="checkbox"/>	Partial completion of one			
			# reps/30 secs		# reps/30 secs				
3	REACH DOWN (reps)	<input type="checkbox"/>	Cannot do at all	Partial = 50% of distance from knee to foot.	<input type="checkbox"/>	Cannot do at all			
		<input type="checkbox"/>	Partial completion of one		<input type="checkbox"/>	Partial completion of one			
			# reps/30 secs			# reps/30 secs			
4	LIFT UP (wt)	Mark highest weight attained:							
		<input type="checkbox"/>	Cannot complete	<input type="checkbox"/>	Lifts ½ kg				
		<input type="checkbox"/>	Lifts empty wrap	<input type="checkbox"/>	Lifts 1 kg				
				<input type="checkbox"/>	Lifts 2 kg				
5	PUSH DOWN (sec)	No partial completion. Max time 30 sec. Record time to nearest 0.1 sec		<input type="checkbox"/>	Cannot complete				
					Seconds held				
6	WRIST UP (reps)	° passive extension		Record degrees of extension where horizontal = 0 Partial = wrist at least to neutral	° passive extension				
		<input type="checkbox"/>	cannot complete		<input type="checkbox"/>	cannot complete			
		<input type="checkbox"/>	Partial completion of one		<input type="checkbox"/>	Partial completion of one			
			Reps/30secs			Reps/30secs			
7a	GRASP DYNAMOMETER	Acquire: yes = able/no = unable to get fingers around dynamometer							
		Release: yes = able/no = unable to extend fingers to let go of dynamometer (not pull arm away)							
		RIGHT			LEFT				
		Acquire		Release		Acquire		Release	
	no	yes	no	yes	no	yes	no	yes	
8	GRASP DYNAMOMETER	RIGHT				LEFT			
		Trial 1	kg:			Trial 1	kg:		
		Trial 2	kg:			Trial 2	kg:		
		Trial 3	kg:			Trial 3	kg:		
9a	LAT. PINCH (credit card) P/F	<input type="checkbox"/>	unable to hold		<input type="checkbox"/>	unable to hold			
		<input type="checkbox"/>	holds for < 5 secs		<input type="checkbox"/>	holds for < 5 secs			
		<input type="checkbox"/>	holds for 5 secs		<input type="checkbox"/>	holds for 5 secs			
9b	LAT. PINCH (force)	Trial 1	kg:		Trial 1	kg:			
		Trial 2	kg:		Trial 2	kg:			
		Trial 3	kg:		Trial 3	kg:			

Subject ID: _____

Item	Item Name	(fill the section below that applies)							
		Right		Both	Left				
10	PULL (wt)	Mark highest weight attained:			Mark highest weight attained:				
		<input type="checkbox"/>	Cannot complete		<input type="checkbox"/>	Cannot complete			
		<input type="checkbox"/>	moves ½ kg		<input type="checkbox"/>	moves ½ kg			
		<input type="checkbox"/>	moves 1 kg		<input type="checkbox"/>	moves 1 kg			
		<input type="checkbox"/>	moves 2 kg		<input type="checkbox"/>	moves 2 kg			
		<input type="checkbox"/>	moves 4 kg	<input type="checkbox"/>	moves 4 kg				
11	PUSH (wt)	Mark highest weight attained:			Mark highest weight attained:				
		<input type="checkbox"/>	Cannot complete		<input type="checkbox"/>	Cannot complete			
		<input type="checkbox"/>	moves ½ kg		<input type="checkbox"/>	moves ½ kg			
		<input type="checkbox"/>	moves 1 kg		<input type="checkbox"/>	moves 1 kg			
		<input type="checkbox"/>	moves 2 kg		<input type="checkbox"/>	moves 2 kg			
		<input type="checkbox"/>	moves 4 kg	<input type="checkbox"/>	moves 4 kg				
7B	Acquire-release part 2 CONTAINER	Acquire: yes = able to get fingers around top of lid; no = unable to get fingers into position							
		Release: yes = able to relax/extend fingers to release lid (not pull arm away); no = unable to extend fingers to release lid							
		RIGHT				LEFT			
		Acquire		Release	Acquire		Release		
		no	yes	no	yes	no	yes		
12	CONTAINER (wt)	RIGHT				LEFT			
		<input type="checkbox"/>	unable to lift			<input type="checkbox"/>	unable to lift		
		<input type="checkbox"/>	lifts empty			<input type="checkbox"/>	lifts empty		
		<input type="checkbox"/>	lifts ½ kg < 5 sec			<input type="checkbox"/>	lifts ½ kg < 5 sec		
		<input type="checkbox"/>	lifts ½ kg for 5 sec			<input type="checkbox"/>	lifts ½ kg for 5 sec		
	<input type="checkbox"/>	lifts 1 kg for 5 sec			<input type="checkbox"/>	lifts 1 kg for 5 sec			
13	PINCH DIE (reps)	<input type="checkbox"/>	Cannot do at all		*must hold with pincer grasp tips of thumb and index finger One rep = both dice moved to other marker and back Partial = both dice moved to other marker	<input type="checkbox"/>	Cannot do at all		
		<input type="checkbox"/>	Partial completion of one			<input type="checkbox"/>	Partial completion of one		
			# reps/30 secs			<input type="checkbox"/>	# reps/30 secs		
14	PENCIL (reps)	<input type="checkbox"/>	Cannot do at all		*must hold with 3-jaw chuck/tripod grasp - thumb, index, middle One rep = pencil to middle slot and back to side slot Partial = pencil to middle slot only	<input type="checkbox"/>	Cannot do at all		
		<input type="checkbox"/>	Partial completion of one			<input type="checkbox"/>	Partial completion of one		
			# reps/30 secs			<input type="checkbox"/>	# reps/30 secs		
15	MANIPULATE (reps)	<input type="checkbox"/>	Cannot do at all		One rep = full 360 degree turn Partial = one half turn	<input type="checkbox"/>	Cannot do at all		
		<input type="checkbox"/>	Partial completion of one			<input type="checkbox"/>	Partial completion of one		
			# reps/30 secs			<input type="checkbox"/>	# reps/30 secs		
16	PUSH INDEX (calculator - sec)	<input type="checkbox"/>	cannot complete		*must use extended index finger Record time to nearest 0.1 second	<input type="checkbox"/>	cannot complete		
			secs to complete			<input type="checkbox"/>	secs to complete		
17	PUSH THUMB (cell phone -sec)	<input type="checkbox"/>	cannot complete		*Must use thumb IP flexion Record time to nearest 0.1 second	<input type="checkbox"/>	cannot complete		
			secs to complete			<input type="checkbox"/>	secs to complete		